

Database, 2018, 1–8 doi: 10.1093/database/bay032 Database update



Database update

Biopanning data bank 2018: hugging next generation phage display

Bifang He^{1,†}, Lixu Jiang^{1,†}, Yaocong Duan¹, Guoshi Chai¹, Yewei Fang¹, Juanjuan Kang¹, Min Yu¹, Ning Li¹, Zhongjie Tang¹, Pengcheng Yao¹, Pengcheng Wu¹, Ratmir Derda² and Jian Huang^{1,*}

¹Center for Informational Biology, University of Electronic Science and Technology of China, No. 2006, Xiyuan Ave, West Hi-Tech Zone, Chengdu 611731, China and ²Department of Chemistry, University of Alberta, 11227 Saskatchewan Drive, Edmonton, AB T6G 2G2, Canada

*Corresponding author: Tel: +86-28-8320-2351; Fax: +86-28-8320-8238; Email: hj@uestc.edu.cn

[†]These authors contributed equally to this work.

Citation details: He,B., Jiang,L., Duan,Y. *et al.* Biopanning data bank 2018: hugging next generation phage display. *Database* (2018) Vol. 2018: article ID bay032; doi:10.1093/database/bay032

Received 8 January 2018; Revised 7 February 2018; Accepted 7 March 2018

Abstract

The 2018 update of the biopanning data bank (BDB) stores phage display data sequenced by Sanger sequencing and next generation sequencing technologies. In this work, we upgraded the database with more biopanning data sets and several new features, including (i) incorporation of next generation biopanning data and the unselected population where the target is not determined and the round of screening is zero; (ii) addition of sequencing information; (iii) improvement of browsing and searching systems and 3 D chemical structure viewer; (iv) integration of standalone tools for target-unrelated peptides analysis within conventional phage display and next generation phage display (NGPD) data. In the current version of BDB (released on 19 January 2018), the database houses 3291 sets of biopanning data collected from 1540 published articles, including 95 NGPD data sets and 3196 traditional biopanning data sets. The BDB database serves as an important and comprehensive resource for developing peptide ligands.

Database URL: The BDB database is available at http://immunet.cn/bdb

Introduction

Phage display allows the high throughput screening of phage-displayed libraries against multiple target molecules such as miRNAs (1), proteins (2, 3), polysaccharides (4), cells (5) and tissues (6). Typically, libraries with billions of peptides are subjected to iterative rounds of affinity selection, commonly referred to as biopanning (7), which rapidly enriches phage clones binding to the target of interest. Due to its high efficiency and versatility, screening of phage-displayed random peptide libraries has found wide applications in the development of diagnostics and therapeutics (8), drug-delivery reagents (9), biomaterials and inorganic functional materials (10, 11).

However, the identification of potential candidates after the biopanning process is a major challenge associated with the phage display technique. Traditionally, 20–100

 $\ensuremath{\mathbb{C}}$ The Author(s) 2018. Published by Oxford University Press.

This is an Open Access article distributed under the terms of the Creative Commons Attribution License (http://creativecommons.org/licenses/by/4.0/), which permits unrestricted reuse, distribution, and reproduction in any medium, provided the original work is properly cited.

individual clones from the final enriched population are randomly picked out and amplified. Then a target-binding screening assay can be used to distinguish the specific binders. Finally, the sequence of the displayed peptide with specific binding could be determined by Sanger sequencing. Given the vast diversity of phage-displayed peptide libraries ($\sim 10^9$ unique sequences), Sanger sequencing analysis is low-throughput and provides a limited perspective (<0.01%) of the complete sequence space. More importantly, traditional phage display is badly hindered by the identification of notorious target-unrelated peptides (TUPs) (12, 13), which may dominate the small sample space and preclude the isolation of target-binding sequences. These TUPs can be divided into two types: selection-related (SrTUPs) and propagation-related (PrTUPs) (14). The selection of SrTUPs is caused by binding to non-target molecules in the biopanning system (15), whereas the isolation of PrTUPs is due to the advantage possessed by certain fast-propagating clones intrinsic to the library (12, 14, 16-18).

To overcome these challenges, many researchers have employed the phage display technology coupled with next generation sequencing (NGPD) in ligand identification (19-24), which allows for a comprehensive characterization of the sequence space of phage-displayed libraries. Next generation sequencing (NGS) technology has also been integrated with other in vitro selection techniques, such as messenger RNA (mRNA) display (25), yeast display (26) and aptamer selection (27). However, next generation phage display (NGPD) also possesses the challenge of evaluating which sequences are the most ideal ligands since the biopanning results are still a mixture of target-binding peptides and TUPs. To exclude any potential TUPs and identify peptides specifically binding to KS483 cells, 'tHoen and coworkers filtered all peptides with a frequency of 2 or higher in the naïve library and hits in PepBank and SAROTUP and selected target-related peptides from the remaining list of peptides (24). They also demonstrated that NGPD can promote the finding of specific binders by alleviating the problem associated with TUPs. Derda et al. have identified glycopeptide ligands for carbohydratebinding proteins by comparing the frequency of each peptide after selection against the target molecule with the corresponding frequency for non-target molecules (28). As can be seen earlier, the strategy to determine target-related peptides can vary from researcher to researcher.

The popularity of the phage display technology inevitably led to the development of databases for phage display data. Artificially selected proteins/peptides database was the first database special for biopanning data (29), which incorporated biopanning data from 195 screening experiments. However, the database has not been upgraded since

2002. To gather these valuable data together, our group implemented and described the MimoDB database (MimoDB 1.0) for peptides selected from phage-displayed random peptide libraries (30). The updated MimoDB 2.0 includes biopanning data isolated from random peptide libraries constructed by diverse display technologies, including phage display, bacterial display, yeast display, ribosome display and mRNA display (31). The database, MimoDB, is actually the abbreviation of mimotope database. As the output of biopanning experiments contains both mimotopes and TUPs, the database was renamed as biopanning data bank (BDB) (32). There are other databases, such as PepBank (33), IEDB (34) and TumorHoPe (35), each harboring a part of biopanning data. These archives are not dedicated to biopanning data storage and take each unique peptide as an individual entry, whereas the result of biopanning usually consists of a set of peptides with different frequencies and affinities. However, it is a remarkable fact that all existing databases are exclusive of NGPD data.

Simultaneously, computational tools for analyzing biopanning data are continuously emerging. These programs can be employed to clean TUPs from the biopanning data (36–38) and interpret biopanning data (39, 40). Unfortunately, all these methods are only applicable to analysis of small-scale biopanning data obtained by traditional phage display. In recent years, powerful analytical methods have been designed for processing and translating sequences, including MATLAB-based single-end and paired-end conversion pipelines (41, 42). Furthermore, software for cluster or motif analysis within massive data has also been proposed, such as multiple specificity identifier (43), Hammock (44) and FASTAptamer (45). Additional methods have been built for finding potential target-binding ligands, such as Enrich2 for any countingbased enrichment/depletion Experiment (46) and PHASTpep for discovery of cell-selective peptides (47). Although these methods are successful in analyzing deep sequencing data, none of them can be utilized to perform TUP analysis within big biopanning data.

In this report, we describe a significant upgrade to the BDB database. We curated NGPD data from peerreviewed publications and integrated standalone data analysis tools, which can facilitate TUP identification within traditional small scale biopanning data and 'big biopanning data' locally. The interactive structure viewer for target-template complex (TTC) or target-peptide complex (TPC) is re-implemented with JSmol and PHP. The search system has also been improved to be more user-friendly and intuitive. With traditional phage display and NGPG data, the BDB database will serve as an important and informational portal for biopanning data and a powerful evidence-based platform for the biopanning community to cross-check their panning results and ensure target-binding specificities. Thus, it provides a comprehensive resource for ligand identification-related studies.

Data upgrade

Inclusion of NGPD and unselected data

As an updated version of BDB released in July 2015 (32) and MimoDB 2.0 (31), the current version of BDB complies with previous literature search methods, data collection approaches, data inclusion criteria and data organization style, which were described in great detail in the literature (31, 32). NGPD screens are more productive than traditional phage display selections, and NGPD data usually have millions of peptide sequences, whereas traditional phage display selections produce several hundreds of peptides. This presents a great challenge as well as an opportunity. In the current version of BDB (released on 19 January 2018), we curated 95 sets of NGPD data from 20 published papers. Sequencing data of selected populations and unselected libraries, for instance, the naïve library and the naïve library after one round or several rounds of amplification were collected. For the unselected dataset, the corresponding target name is not determined and the round of panning is zero (see Figure 1).

ngDataSet Target Template Library Structu

New data fields

With the advance of sequencing technology, it is necessary to include sequencing information into BDB. For each entry of the biopanningdataset table, additional fields, including sequencing platform and sequencing method, were curated (see Figure 1, indicated by red boxes). According to our records, the BDB database stores peptide data sequenced by various high-throughput sequencing platforms, including Illumina Genome Analyzer[®], Roche 454 GS System[®], Ion Torrent system[®] and Applied Biosystems ABI PRISM system®. However, most published articles using Sanger sequencing did not provide any information about sequencing platforms. Therefore, for most Sanger sequencing data sets, corresponding sequencing systems are unknown. We also added curation information, including curator and curation date (see Figure 1, marked by red boxes).

Data summary and statistics

The BDB database released recently (19 January 2018) contains a total of 3291 sets of phage display data manually curated from 1540 peer-reviewed articles, a significant increase over the 2904 biopanning data sets in the previous major release (22 July 2015) (32). According to the sequencing technology, these data sets are distributed into 3196 sets of traditional biopanning data and 95 sets of NGPD data.

ngDataSet ID: Peptides: WALATW (1 VNAATC (12 LCAEAC (10 nces Features Motif: 20 imental Method encing Platform: Illumina HiSeq 2000 and/or Genome Analyzer Iix Next-generation sequencing Method: nds of Pa PMID: 22870226 Not determined Not determined Not determined Template Name: Structure of Target-Template Complex: Not determined Structure of Target-Peptide Complex: Not determined Library Name: Affinity Measurement Method: Affinity Measurement Descript pVIII fth1-dp-based X7 phage display library (X7) ent Descriptio For the analysis of the naive libraries, four samples of 4.0e8 phages each were amplified individually by the PCR reactions. The amplified PCR products were validated for size (152 bp) by running in 2% agarose gel. All 4 PCR cleaned products were united and their concentration was measured. The sample was dired by speed vas: and sert for illumina sequencing. 4 PCR cleaned products were united and their concentration was measured. The The 20 top most frequent peptides for each capture are given. Brief Descriptio ningDataSet Comments: Curator Lixu Jiang **Curation Date** 2017-05-09

Figure 1. An example of a set of phage display data from the naïve library. New data fields are indicated by red boxes. For the unselected dataset, the corresponding target name is not determined and the round of panning is zero (highlighted by blue boxes).

Based on their target types, biopanning data sets can be divided into nine categories. As shown in Figure 2A, we can conclude that proteins, cells and inorganic molecules or substances are the three most commonly used targets in phage display screening experiments. Importantly, the current release of BDB also contains 29 795 peptides, 1971 targets, 545 templates, 484 different combinatorial libraries and 318 three-dimensional structures of TTC or TPC resolved experimentally (for instance, by X-ray crystallography or nuclear magnetic resonance spectroscopy). Generally, the BDB database is revised and updated once per quarter. Since its first release on 30 August 2010, the database has been updated as many as 29 times. We summarized the data entries increase annually (Figure 2B).

Utilities upgrade

Browsing BDB

As described previously, each table in the BDB database can be quickly browsed via the 'Browse' drop-down menu or the 'Browse' option on the left side of the search panel (32). To improve user experience, we redesigned the BDB homepage layout. Currently, traditional or NGPD data can be browsed directly by clicking the corresponding figure on the left side of the BDB homepage (Figure 3). Thus users can conveniently visit these two types of data separately. When browsing the biopanningdataset table, if the entry contains >200 peptide sequences, only the first 10 peptides will be displayed and the whole data set can be downloaded via the provided link (Figure 4). It is particularly worth mentioning that the on-the-fly data visualization tool, sequences features, is disabled when a set of phage display data contains <3 or >200 peptides (Figure 4).

Searching BDB

Based on user feedback, more data fields, such as sequencing platform, sequencing method, number of unique peptides, have been added to the advanced searching system. In the current release of BDB, most data fields in the database can be searched through the advanced searching system. For each data field to be searched, the search program will provide prompt and appropriate hint for users and help them quickly find data of interest, leading to an improved user experience.

Viewing complex structures

Java Applets have served as a primary program in displaying interactive 3 D chemical structures on a web page for about two decades. In fact, Jmol (Jmol: an open-source Java viewer for chemical structures in 3 D, http://www. jmol.org/) has been the 3 D viewer for the BDB website



Figure 2. Summary of data entries in the BDB database. (A) Biopanning data sets are grouped into nine classes according to their target types. The number of biopanning data set selected by each type of target is shown. (B) Annual increase of data entries in BDB.



Figure 3. A screenshot of the BDB homepage. The homepage of BDB clearly shows that the BDB database contains both traditional and NGPD data. Users can browse these two kinds of data via clicking the corresponding figure on the left side of the website.

Biopanning Data	
biopanning bata	
BiopanningDataSet ID: 2210	
Peptides: Sequences	
VFIPYGH	
KWAGTHI	
NV KALANI NU KALANI	
HAP LOON	
RDLATLO	
TVINENY	
GWWHLPR	
Only first 10 sequences were listed above. For more data, please Download, tgz file.	
Sequences Features:	
Motif	
Unique Sequence Number: 254	
Experimental Method: Phage display (competitive panning)	
Sequencing Platform: Jumina Hiseq 2000	
Sequencing Method: Next-generation sequencing	
Rounds of Panining: 2 Peferanatic Paline 2006227	
Tarnot Name: Anti-DAD cora cample / DAD/	
Tamplate Name: Destatic acid historiataca DAD	
Structure of Tarret-Template Complex: Not determined	
Structure of Target-Pentide Complex: Not determined.	
Library Name: Ph.D7 phage display library (X7)	
Affinity Measurement Method:	
Affinity Measurement Description:	
Brief Description: The phage library was mixed with mouse or human serum. After incubation and washing, protein G was added to the phage-antibody mixture.	
BiopanningDataSet Comments:	
Curator:	
Curation Date: 0000-00-00	

Figure 4. An example of a set of biopannning data with >200 peptides. If a set of biopanning data contains >200 peptides, the sequences features tool is disabled, and only the first 10 peptide sequences will be displayed and the complete data set is available for download via the provided link.

(30–32) and other web resources, such as the research (14) for structural bioinformatics protein data bank (RCSB PDB, http://rcsb.org) (48). Nevertheless, several Internet browsers (e.g. Chrome) no longer support Java Applets because of concerns for security precautions. In 2016, Oracle announced that Java Applets would be deprecated with Java version 9.0. To provide improved support for 3 D

graphics, the BDB website has implemented JSmol (JSmol: an open-source HTML5 viewer for chemical structures in 3 D, http://wiki.jmol.org/index.php/JSmol) to interactively view TTC or TPC, which enables Jmol to display 3 D molecular structures on devices that do not have Java installed, or for which Java is not available or has not been installed.

Using BDB-powered tools

In last major upgrade of BDB, we integrated a panel of tools in the SAROTUP suite into BDB (32). These tools can be classified into database-based, motif-based and machine learning methods-based tools and have played a significant role in cleaning TUPs from phage-displayed libraries and one-bead-one-compound combinatorial peptide libraries (49-52). However, these tools are unable to analyze next generation sequencing biopanning data. In the current release of BDB, we integrated standalone graphical user interface (GUI) and command-line versions for each tool developed using open source Qt 5.6 and C++ language, which can be used to find TUPs within biopanning results of conventional phage display and NGPD, whereas the web-based tools can be only employed to handle small-scale data. The offline version of each tool can run locally in Windows or Linux systems. The interface and utilization of the GUI version are similar to that of the web server. We also added a new machine learning method-based tool, known as PSBinder (53). This program allows users to predict putative polystyrene surfacebinding peptides from biopanning data or to find novel candidates for polystyrene affinity tags. All tools can be available for download and free use. As shown in Figure 5A, six tools, i.e. TUPScan, MimoSearch,

MimoBlast, PhD7Faster, SABinder and PSBinder, can be employed to perform TUP analysis, while MimoScan allows users to check how specific their consensus sequences, motifs or patterns derived from panning results are. We take MimoSearch as an example to show the utilization of these tools (Figure 5B). All these tools are valuable to the biopanning community. With the accumulation of more biopanning data in BDB, the database-based tools, including MimoSearch, MimoBlast, MimoScan, will be increasingly powerful and popular. We expect that cleaning TUPs from biopanning data will be adopted as a necessary and standard procedure in ligand development.

Conclusion and future development

The BDB database is the largest archive for biopanning data, which provides a valuable resource for ligand development related studies. Currently, the database incorporates phage display data sequenced by both Sanger sequencing and high-throughput sequencing platforms. With the increasing popularity of NGPD, more sets of NGPD data will be added into the BDB database. The database will be continually maintained and updated. Furthermore, an experimentally supported dataset will be constructed and available in the near future, which will



Figure 5. BDB-powered tools. (A) Seven tools surrounded by a red box, i.e. MimoSearch, MimoBlast, MimoScan, TUPScan, PhD7Faster, SABinder and PSBinder, can be accessible by clicking the secondary menu items from the 'Tools' drop-down menu. (B) We take MimoSearch as an example to show the utilization of these tools. Users can enter a set of peptide sequences in the text area or upload a sequence file. After submission, the result will be returned and displayed in a table.

provide candidate peptides for drug and vaccine design and therapeutics and diagnostics development. Computational tools for clustering analysis, motif analysis and epitope mapping will be integrated into the database.

Acknowledgements

The authors are grateful to the anonymous reviewers for their valuable suggestions and comments, which have led to the improvement of this paper.

Funding

This work was supported by the National Natural Science Foundation of China [grant no. 61571095] and the Fundamental Research Funds for the Central Universities of China [grant no. ZYGX2015Z006]. Funding for open access charge: National Natural Science Foundation of China [grant no. 61571095].

Conflict of interest. None declared.

References

- Bose, D., Nahar, S., Rai, M.K. *et al.* (2015) Selective inhibition of miR-21 by phage display screened peptide. *Nucleic Acids Res.*, 43, 4342–4352.
- He,B., Mao,C., Ru,B. *et al.* (2013) Epitope mapping of metuximab on CD147 using phage display and molecular docking. *Comput. Math. Methods Med.*, 2013, 1.
- Zhang,Y., He,B., Liu,K. *et al.* (2017) A novel peptide specifically binding to VEGF receptor suppresses angiogenesis *in vitro* and *in vivo*. *Signal Transduct. Target. Ther.*, 2, 17010.
- 4. Agrawal, S., Kulabhusan, P.K., Joshi, M. *et al.* (2016) A high affinity phage-displayed peptide as a recognition probe for the detection of *Salmonella typhimurium*. *J Biotechnol.*, **231**, 40–45.
- Gray, B.P. and Brown, K.C. (2014) Combinatorial peptide libraries: mining for cell-binding peptides. *Chem. Rev.*, 114, 1020–1081.
- Yang,Y.H., Rajaiah,R., Ruoslahti,E. *et al.* (2011) Peptides targeting inflamed synovial vasculature attenuate autoimmune arthritis. *Proc Natl Acad Sci USA*, 108, 12857–12862.
- Pande, J., Szewczyk, M.M. and Grover, A.K. (2010) Phage display: concept, innovations, applications and future. *Biotechnol. Adv.*, 28, 849–858.
- Gautam,A., Kapoor,P., Chaudhary,K. *et al.* (2014) Tumor homing peptides as molecular probes for cancer therapeutics, diagnostics and theranostics. *Curr. Med. Chem.*, 21, 2367–2391.
- Liu, R., Li, X., Xiao, W. et al. (2016) Tumor-targeting peptides from combinatorial libraries. Adv. Drug. Deliv. Rev., 110–111, 13–37.
- Mannocci,L., Leimbacher,M., Wichert,M. et al. (2011) 20 years of DNA-encoded chemical libraries. Chem. Commun. (Camb), 47, 12747–12753.
- 11. Lee, Y.J., Yi, H., Kim, W.J. *et al.* (2009) Fabricating genetically engineered high-power lithium-ion batteries using multiple virus genes. *Science*, **324**, 1051–1055.
- Matochko, W.L., Cory Li, S., Tang, S.K. *et al.* (2014) Prospective identification of parasitic sequences in phage display screens. *Nucleic Acids Res.*, 42, 1784–1798.

- Vodnik, M., Zager, U., Strukelj, B. *et al.* (2011) Phage display: selecting straws instead of a needle from a haystack. *Molecules*, 16, 790–817.
- 14. Thomas, W.D., Golomb, M. and Smith, G.P. (2010) Corruption of phage display libraries by target-unrelated clones: diagnosis and countermeasures. *Anal. Biochem.*, 407, 237–240.
- Vodnik, M., Strukelj, B. and Lunder, M. (2012) HWGMWSY, an unanticipated polystyrene binding peptide from random phage display libraries. *Anal. Biochem.*, 424, 83–86.
- Nguyen,K.T., Adamkiewicz,M.A., Hebert,L.E. *et al.* (2014) Identification and characterization of mutant clones with enhanced propagation rates from phage-displayed peptide libraries. *Anal. Biochem.*, 462, 35–43.
- Zygiel,E.M., Noren,K.A., Adamkiewicz,M.A. *et al.* (2017) Various mutations compensate for a deleterious lacZalpha insert in the replication enhancer of M13 bacteriophage. *PLoS One*, 12, e0176421.
- Brammer,L.A., Bolduc,B., Kass,J.L. *et al.* (2008) A target-unrelated peptide in an M13 phage display library traced to an advantageous mutation in the gene II ribosome-binding site. *Anal. Biochem.*, 373, 88–98.
- Rentero Rebollo,I., Sabisz,M., Baeriswyl,V. *et al.* (2014) Identification of target-binding peptide motifs by highthroughput sequencing of phage-selected peptides. *Nucleic Acids Res.*, 42, e169.
- 20. Ryvkin, A., Ashkenazy, H., Smelyanski, L. *et al.* (2012) Deep panning: steps towards probing the IgOme. *PLoS One*, 7, e41469.
- Ngubane, N.A., Gresh, L., Ioerger, T.R. *et al.* (2013) High-throughput sequencing enhanced phage display identifies peptides that bind mycobacteria. *PLoS One*, 8, e77844.
- Ravn,U., Didelot,G., Venet,S. *et al.* (2013) Deep sequencing of phage display libraries to support antibody discovery. *Methods*, 60, 99–110.
- Staquicini,F.I., Cardo-Vila,M., Kolonin,M.G. *et al.* (2011) Vascular ligand-receptor mapping by direct combinatorial selection in cancer patients. *Proc. Natl. Acad. Sci. USA*, 108, 18637–18642.
- 24. t Hoen,P.A., Jirka,S.M., Ten Broeke,B.R. *et al.* (2012) Phage display screening without repetitious selection rounds. *Anal. Biochem.*, 421, 622–631.
- Olson, C.A., Nie, J., Diep, J. *et al.* (2012) Single-round, multiplexed antibody mimetic design through mRNA display. *Angew. Chem. Int. Ed. Engl.*, 51, 12449–12453.
- Van Blarcom, T., Rossi, A., Foletti, D. *et al.* (2015) Precise and efficient antibody epitope determination through library design, yeast display and next-generation sequencing. *J. Mol. Biol.*, 427, 1513–1534.
- 27. Schutze, T., Wilhelm, B., Greiner, N. *et al.* (2011) Probing the SELEX process with next-generation sequencing. *PLoS One*, *6*, e29604.
- Ng,S., Lin,E., Kitov,P.I. *et al.* (2015) Genetically encoded fragment-based discovery of glycopeptide ligands for carbohydrate-binding proteins. *J. Am. Chem. Soc.*, 137, 5248–5251.
- Valuev,V.P., Afonnikov,D.A., Ponomarenko,M.P. *et al.* (2002) ASPD (artificially selected proteins/peptides database): a database of proteins and peptides evolved in vitro. *Nucleic Acids Res.*, 30, 200–202.

- Page 8 of 8
- Ru,B., Huang,J., Dai,P. *et al.* (2010) MimoDB: a new repository for mimotope data derived from phage display technology. *Molecules*, 15, 8279–8288.
- 31. Huang, J., Ru, B., Zhu, P. et al. (2012) MimoDB 2.0: a mimotope database and beyond. Nucleic Acids Res., 40, D271–D277.
- 32. He,B., Chai,G., Duan,Y. *et al.* (2016) BDB: biopanning data bank. *Nucleic Acids Res.*, 44, D1127–D1132.
- Shtatland, T., Guettler, D., Kossodo, M. *et al.* (2007) PepBank—a database of peptides based on sequence text mining and public peptide data sources. *BMC Bioinformatics*, 8, 280.
- Fleri, W., Paul, S., Dhanda, S.K. *et al.* (2017) The immune epitope database and analysis resource in epitope discovery and synthetic vaccine design. *Front. Immunol.*, 8, 278.
- 35. Kapoor, P., Singh, H., Gautam, A. *et al.* (2012) TumorHoPe: a database of tumor homing peptides. *PLoS One*, 7, e35187.
- Huang, J., Ru, B., Li, S. *et al.* (2010) SAROTUP: scanner and reporter of target-unrelated peptides. *J. Biomed. Biotechnol.*, 2010, 1.
- Ru,B., T Hoen,P.A., Nie,F. *et al.* (2014) PhD7Faster: predicting clones propagating faster from the Ph.D.-7 phage display peptide library. *J. Bioinform. Comput. Biol.*, 12, 1450005.
- He,B., Kang,J., Ru,B. *et al.* (2016) SABinder: a web service for predicting streptavidin-binding peptides. *Biomed. Res. Int.*, 2016, 9175143.
- Negi,S.S. and Braun,W. (2009) Automated detection of conformational epitopes using phage display peptide sequences. *Bioinform. Biol. Insights*, 3, BBI.S2745–B81.
- Huang, J., Gutteridge, A., Honda, W. et al. (2006) MIMOX: a web tool for phage display based epitope mapping. BMC Bioinformatics, 7, 451.
- Matochko, W.L., Chu, K., Jin, B. *et al.* (2012) Deep sequencing analysis of phage libraries using Illumina platform. *Methods*, 58, 47–55.
- 42. He,B., Tjhung,K.F., Bennett,N.J. *et al.* (2018) Compositional bias in naive and chemically-modified phage-displayed libraries uncovered by paired-end deep sequencing. *Sci. Rep.*, 8, 1214.
- 43. Kim, T., Tyndel, M.S., Huang, H. et al. (2012) MUSI: an integrated system for identifying multiple specificity from very

large peptide or nucleic acid data sets. Nucleic Acids Res., 40, e47.

- 44. Krejci, A., Hupp, T.R., Lexa, M. et al. (2016) Hammock: a hidden Markov model-based peptide clustering algorithm to identify protein-interaction consensus motifs in large datasets. *Bioinformatics*, 32, 9–16.
- Alam,K.K., Chang,J.L. and Burke,D.H. (2015) FASTAptamer: a bioinformatic toolkit for high-throughput sequence analysis of combinatorial selections. *Mol. Ther. Nucl. Acids*, 4, e230.
- Rubin,A.F., Gelman,H., Lucas,N. *et al.* (2017) A statistical framework for analyzing deep mutational scanning data. *Genome Biol.*, 18, 150.
- Brinton,L.T., Bauknight,D.K., Dasa,S.S. *et al.* (2016) PHASTpep: analysis software for discovery of cell-selective peptides via phage display and next-generation sequencing. *PLoS One*, **11**, e0155244.
- Rose, P.W., Prlic, A., Altunkaya, A. *et al.* (2017) The RCSB protein data bank: integrative view of protein, gene and 3D structural information. *Nucleic Acids Res.*, 45, D271–D281.
- Hou,L., Zhu,D., Liang,Y. *et al.* (2017) Identification of a specific peptide binding to colon cancer cells from a phage-displayed peptide library. *Br. J. Cancer*, 118, 79–87.
- Christiansen, A., Kringelum, J.V., Hansen, C.S. *et al.* (2015) High-throughput sequencing enhanced phage display enables the identification of patient-specific epitope motifs in serum. *Sci. Rep.*, 5, 12913.
- Liu,G.W., Livesay,B.R., Kacherovsky,N.A. *et al.* (2015) Efficient identification of murine M2 macrophage peptide targeting ligands by phage display and next-generation sequencing. *Bioconjug. Chem.*, 26, 1811–1817.
- Leung,N.Y., Wai,C.Y., Ho,M.H. *et al.* (2017) Screening and identification of mimotopes of the major shrimp allergen tropomyosin using one-bead-one-compound peptide libraries. *Cell Mol. Immunol.*, 14, 308–318.
- Li,N., Kang,J., Jiang,L. *et al.* (2017) PSBinder: a web service for predicting polystyrene surface-binding peptides. *Biomed. Res. Int.*, 2017, 1.